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REFERENTIES

Adcock DM, Kressin DC, Marlar RA.

Effect of 3.2% vs 3.8% sodium citrate concentration on routine coagulation testing. *Am J Clin Pathol* 1997; 107: 105-110.

Adcock DM, Kressin DC, Marlar RA.

The effect of time and temperature variables on routine coagulation tests. *Blood Coagul Fibrinolysis* 1998; 9: 463-470.

Adcock DM, Kressin DC, Marlar RA.

Minimum specimen volume requirements for routine coagulation testing. Dependence on citrate concentration. *Am J Clin Pathol* 1998; 109:595-599.

Adcock Funk DM, Lippi G, Favaloro EJ.

Quality standards for sample processing, transportation, and storage in hemostasis testing. *Semin Thromb Hemost* 2012; 38:1-10.

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Brien WF, Schaus MR, Cooper KE, O'Keefe BT, Inwoon M.

Lupus anticoagulant testing: effect of platelet count on the activated partial thromboplastin time. *Br J Biomed Sc* 1993; 50: 114-116.

Brandt JT, Triplett DA, Alving B, Scharrer I.

Criteria for the diagnosis of lupus anticoagulants: an update. On behalf of the Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibody of the Scientific and Standardisation Committee of the ISTH. *Thromb Haemost* 1995; 74:1185-1190.

Capel P, Chatelain B, Leclercq R, Lust A, Masure R, Arnout J.

Quality control in haemostasis. *Acta Clinica Belgica* 1992; 47: 308-318.

CLSI: Procedures for the collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-sixth Edition. *CLSI document H3-A6.2007*

CLSI guideline:

Collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays; Approved guideline – Fifth Edition. H21-A5. Wayne, PA, USA: Clinical and Laboratory Standards Institute; 2008.

Contant G, Gouault-Heilmann M, Martinoli JL.

Heparin inactivation during blood storage: its prevention by blood collection in citric acid, theophylline, adenosine, dipyridamole –C.T.A.D. mixture-. *Thromb Res* 1983; 31:365-374.

Feng LM, Zhao Y, Zhao HC, Shao ZX.

Effects of storage time and temperature on coagulation tests and factors in fresh plasma. *Sci Rep* 2014; 4:3868

Fiebig EW, Etzell JE, Ng VL.

Clinically relevant differences in prothrombin time and INR values related to blood sample collection in plastic vs glass tubes. *Am J Clin Pathol.* 2005;124:902-909.

Gogstad GO, Dale S, Brosstad F, Brandsnes Ø, Holtlund J, Mørk E, Gärtner E, Borch SM.

Assay of D-dimer based on immunofiltration and staining with gold colloids. *Clin Chem* 1993 ; 39 :2070-2076.

Greaves M, Cohen H, Machin SJ, Mackie I.

Guidelines on the investigation and management of the antiphospholipid syndrome. *Br J Haematol* 2000; 109: 704-715.

Ho CH, Wu S-Y.

The influence of time, temperature and packed cell on activated partial thromboplastin time and prothrombin time. *Thromb Res* 1991; 62: 625-633.

Hope E, Mayorga SR, Robinson S, Goldberg L, Leveson JE, Marengo-Rowe AJ, Peerschke EIB.

Preparation of plasma for coagulation testing. Evaluation of the StatSpin high-speed centrifuge. *Laboratory Medicine* 1991; 22: 190-193.

Hillyard CJ, Blake AS, Wilson K, Rylatt DB, Miles S, Bunch R, Elms MJ, Barnes A, Bundesen PG.

A latex agglutination assay for D dimer: evaluation and application to the diagnosis of thrombotic disease. *Clin Chem* 1987; 33:1837-1840.

Iversen LH.

Pre-analytical variation in the measurements of sensitive markers of coagulation and fibrinolysis: the influence of venipuncture and mixing of blood.

Haemostasis 1997; 27: 119-124.

Kemkes-Matthes B, Fischer R, Peetz D.

Influence of 8 and 24-h storage of whole blood at ambient temperature on prothrombin time, activated partial thromboplastin time, fibrinogen, thrombin time, antithrombin and D-dimer.

Blood Coagul Fibrin 2011;22:215-220.

Laga AC, Cheves TA, Sweeney JD.

The effect of specimen hemolysis on coagulation test results.

Am J Clin Pathol 2006; 126:748-755.

Lippi G, Montagnana M, Salvagno GL, Guidi GC.

Interference of blood cell lysis on routine coagulation testing.

Arch Pathol Lab Med 2006; 130:181-184.

Lippi G, Salvagno GL, Montagnana M, Guidi GC.

Short-term venous stasis influences routine coagulation testing. *Blood Coagul Fibrinolysis* 2005; 16: 453-458.

Lippi G, Salvagno GL, Poli G, Guidi GC.

Influence of centrifugation temperature on routine coagulation testing. *Clin. Chem.* 2006; 52: 537-538.

Molenaar PJ, Leyte A.

Pre-acquisition system assessment of the Sysmex(®) Coagulation System CS-2100i and comparison with end-user verification; a model for the regional introduction of new analysers and methods. *Clin Chem Lab Med*

2011; 49: 1479-89

Peterson P, Gottfried EL.

The effects of inaccurate blood sample volume on prothrombin time (PT) and activated partial thromboplastin time (aPTT). *Thromb Haemost* 1982; 47:101-103.

Pretorius L, Janse van Rensburg WJ, Conradie C, Coetzee MJ.

Minimum citrate tube fill volume for routine coagulation testing.. *Int J Lab Hemat* 2014; 36: 493-495.

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Collection of blood specimens by venipuncture for plasma-based coagulation assays.

Necessity of a discard tube.

Am J Clin Pathol 2010; 133: 331-335

Rao LV, Okorodudu AO, Petersen JR,

Stability of prothrombin time and activated partial thromboplastin time tests under different storage conditions.

Clin Chim Acta 2000; 300: 13-21

Ray MJ, Carroll PA, Just SJE, Hanson GAT.

A low volume specimen container suitable for monitoring the aPTT of heparinized patients. *Blood Coagul Fibrinol* 1993; 4: 805-807.

Reneke J, Ezzell J, Leslie S, Ng VL, Gottfried EL.

Prolonged prothrombin time and activated partial thromboplastin time due to underfilled specimen tubes with 109 mmol/l (3.2%) citrate anticoagulant.

Am J Clin Pathol 1998; 109: 754-757.

Sletnes KE, Gravem K, Wisloff F.

Preparation of plasma for the detection of lupus anticoagulants and antiphospholipid antibodies. *Thromb Res* 1992; 66: 43-53.

Van den Besselaar AMHP, Chantarangkul V, Tripodi A.

A comparison of two sodium citrate concentrations in two evacuated blood collection systems for prothrombin time and ISI determination. *Thromb Haemost* 2000 ; 84 :664-667.

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Effect of magnesium contamination in evacuated blood collection tubes on the prothrombin time test and ISI calibration using recombinant human thromboplastin and different types of coagulometer. *Thromb Res* 2005 ; 115 :239-244.

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Comparative study of blood collection tubes and thromboplastin reagents for correction of INR discrepancies: a proposal for maximum allowable magnesium contamination in sodium citrate anticoagulant solutions.

Am J Clin Pathol 2012; 138:248-254.

Van Geest-Daalderop JHH, Mulder AB, Boonman-de Winter LJM, Hoekstra MMCL, Van den Besselaar AMHP.

Preanalytical variables and off-site blood collection : influences on the results of the prothrombin time/international normalized ratio test and implications for monitoring of oral anticoagulant therapy. *Clin Chem* 2005 ; 51 :561-568.

Ver Elst K, Vermeiren S, Schouwens S, Callebaut V, Thomson W, WeekX S. Validation of minimal citrate tube fill volume routine coagulation test on ACL TOP 500 CTS. *Int Jnl. Lab. Hem.* 2013; 35: 614-619

Walker ID.

Blood collection and sample preparation: pre-analytical variation. In: J Jespersen, RM Bertina, F Haverkate (editors), 2nd revised edition of ECAT assay procedures. *A manual of laboratory techniques. Kluwer Academic Publishers, Dordrecht, 1999; pp 21-28*

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Woodhams B, Girardot O, Blanco MJ, Colesse G, Gourmelin.

Stability of coagulation proteins in frozen plasma. *Blood Coagul Fibrinolysis* 2001; 12:229-236.

Zürcher M, Sulzer I, Barizzi G, Lämmle B, Alberio L.

Stability of coagulation assays performed in plasma from citrated whole blood transported at ambient temperature. *Thromb Haemost* 2008; 99: 416-426.

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
$0,20 - 0,13 = 0,07$ ml citraat per ml plasma.

Reagents

Select the curve for the total volume of anticoagulated blood required (10, 5, 2, or 1 mL). Enter the chart at the patient's packed cell volume on the horizontal axis and read off the corresponding volume of anticoagulant solution on the vertical axis. Place the volume in a collection tube and add blood up to the required total volume. The chart assumes a normal packed cell volume of 43%; in that case, 1 volume of anticoagulant is made up to 10 volumes with blood.

Alternatively, the anticoagulant volume may be calculated from the expression:

$$X = (100 - PCV) / (595 - PCV) \text{ vol,}$$

where, X is the volume of anticoagulant required to prepare unit volume of anticoagulated blood, and PCV is the packed cell volume in %; e.g., to determine the volume required for 5 mL anticoagulated blood, calculate 5 x mL.

(Reproduced from Ingram, G.I.C., Brozovic, M., and Slater, N.G.P. Bleeding disorders, investigations and management, 2nd edition. Blackwell's Scientific Publications, Oxford: pp. 244-245, 1982, by permission of the authors and publishers.)

CLSI H21-A5