



Diagnostic evaluation of the serological detection performance on HBV surface antigens in the Netherlands

Authors:

E. Vanhulle^{1,2}, M. Simões^{1,2}, H. Ten Bruin³, G. Boland⁴, A. Brandenburg⁵, and K. Benschop¹; on behalf of participating laboratories⁶

Affiliations:

¹ Centre for Infectious Disease Research, Diagnostics and Laboratory Surveillance (IDS), RIVM, Bilthoven, the Netherlands.

² European Program for Public Health Microbiology Training (EUPHEM), European Centre for Disease Prevention and Control (ECDC), Stockholm, Sweden.

³ MCA Laboratorium, Streekziekenhuis Koningin Beatrix, Winterswijk, the Netherlands.

⁴ University Medical Center, Utrecht, the Netherlands.

⁵ Certe medische microbiologie Friesland en NOP, Leeuwarden, the Netherlands.

⁶ Contact persons of participating MMLs:

Andeweg K., Reinier Haga MDC; Tjong Ayong N.M., Health Care Laboratory St. Maarten NV; Berkhout H., Dicoon Canisius-Wilhelmina Ziekenhuis; Cremers-van der Crommert K., Eurofins PAMM B.V.; Crielaard A., MML IJsselland Ziekenhuis; Daemen P., Radboudumc; de Vries J.J.C., Leiden University Medical Center; Herpers B., Stichting Streeklaboratorium voor de Volksgezondheid Kennemerland; Heusinkveld M., Ziekenhuis Gelderse Vallei; Hira V., Groene Hart Ziekenhuis; Hoogewerf M., Noordwest Ziekenhuis groep; Huijgen H.J., Rode Kruis Ziekenhuis; Kalpoe J., LABHOH (Aruba); Kissing J., Eurofins Gelre; Koene F., Medical Laboratory Services Curaçao; Meulenberg D., Admiraal de Ruyter Ziekenhuis; Priem I., Saltro B.V.; Raes-van Heerenbeek K., Diagnostiek Voor U; Riezebos-Brilman A., Laboratory for Medical Microbiology and Public Health; Roos A., St. Antonius Ziekenhuis; Scholing M., OLVG Lab B.V.; van Corteveen-Splinter A., CBSL TergooiMC; van der Sluijs L., Haga Ziekenhuis; van Tienen C., MaasstadLab LMM; Velzing J., Erasmus MC Rotterdam; Verduyn Lunel F., UMC Utrecht.

*Contact: kim.benschop@rivm.nl



Contents

Abstract.....	1
1. Introduction	2
HBV genotype F	3
HBV diagnosis	4
HBV surveillance in the Netherlands	4
Role of EQA/SKML.....	6
Objectives	7
2. Study design and methods.....	8
Selection and characteristics of the EQA subgenotypes	8
Procedure for participating laboratories.....	9
Identification of eligible laboratories.....	9
Distribution of EQA subgenotypes to laboratories and reporting of results.....	9
Analysis of the results.....	10
Fictional standardization.....	10
3. Results.....	12
Observed challenges.....	19
4. Conclusions	20
5. Recommendations	21
Acknowledgments.....	23
Disclaimer.....	23
References	24
Appendices.....	26
Appendix 1. Invitation letter.....	26
Appendix 2. Reporting form.	28
Appendix 3. Instructions for use of the trial EQA panel.	33
Appendix 4. Steps followed for fictional standardization.....	34

Supplementary Figure 1..... 35

Supplementary Figure 2..... 36

1 Abstract

2 Hepatitis B virus (HBV) poses a substantial global health challenge, with close to 300 million individuals
3 experiencing chronic infections worldwide. Categorized into ten distinct types (A-J), HBV exhibits a
4 high genetic variability. In the Netherlands, genotype A is most prevalent, followed by genotype D.
5 However, since 2019, a notable increase in genotype F is reshaping the epidemiological landscape of
6 HBV infections. The impact of HBV's high genetic variability on serological assay performance,
7 especially in detection efficacy of HBV surface antigen (HBsAg) from different genotypes by assays
8 used in routine diagnostics, remains an understudied area. Here, we report on the results obtained
9 from the study that aimed to investigate the serological detection performance on HBsAg diagnostic
10 assays used in Dutch Medical Microbiology and Clinical Chemical laboratories for different HBV
11 genotypes.

12 An external quality assessment (EQA) panel containing sixteen well-characterized serological samples
13 with known antigen amounts, representing various HBV genotypes (A, B, C, D, E, F, and H; WHO
14 international HBV reference panel (1)), including a negative serum control, was designed and
15 distributed to twenty-seven participating Dutch laboratories. Each sample contained between 21 and
16 32 IU/mL of HBsAg. Laboratories were instructed to perform routine diagnostic tests and report their
17 results.

18 Three quantitative and eight qualitative diagnostic assays were evaluated. The EQA data revealed
19 consistent detection of all genotypes, including genotype F2. Quantitative assays demonstrated
20 variability in antigen detection of the same genotype among laboratories employing the same
21 diagnostic system, e.g., Liaison XL (DiaSorin). Additionally, all quantitative assays exhibited variations
22 in HBsAg detection rates for different HBsAg genotypes, with values ranging from 10 IU/mL
23 (subgenotype A2) to 76 IU/mL (subgenotype E). The lowest measured values (average HBsAg < 29
24 IU/mL) were observed for genotypes A2, B2, D1, D2, D3, and F2.

25 Our study shows consistent serological detection of HBsAg of the evaluated distinct HBV genotypes,
26 indicating good assay performance of the tested immunoassays. However, variations in detection
27 levels among different genotypes, and assay performance discrepancies have been observed and
28 warrant attention with diagnostic interpretation of HBsAg values, in particular for antiviral resistance
29 testing.

30 1. Introduction

31 According to the latest available data, Hepatitis B virus (HBV) poses a substantial global health
32 challenge, with an estimated 250 - 260 million individuals experiencing chronic infections worldwide.
33 Annually, approximately 1.1 million people die due to complications associated with HBV (2). The
34 WHO has initiated efforts towards elimination of HBV by 2030 by active case finding, contact training,
35 treatment, and vaccination plans.

36 Characterized as a small, enveloped, double-stranded DNA reverse-transcribing virus
37 (*Orthohepadnavirus*, *Hepadnaviridae* family), HBV exhibits a preference for hepatocytes, leading to
38 both acute and chronic hepatic infections. The primary mode of HBV transmission is vertical
39 transmission (3). Other modes of transmission include blood and sexual transmission, as well as
40 intrafamilial spread, which may occur through prolonged contact between children or household
41 members, likely via exposure to fluids from infected individuals (4). HBV is known to exhibit a high
42 genetic variability because of the lack of the proof-reading mechanism (5). Therefore HBV is
43 categorized into 10 distinct types (A-J) with a genome-wide intergenomic sequence divergence of at
44 least 7.5% (6). Additionally, HBV is categorized into four major serotypes (adr, adw, ayr, ayw) based
45 on antigenic epitopes found on its envelope proteins. These serotypes are characterized by a common
46 determinant present in all wild-type variants (a) and two mutually exclusive determinant allele pairs
47 (d/y and w/r) (7). The geographic distribution of HBV genotypes (Table 1) is well-defined and correlates
48 strongly with regional host populations, as well as socio-demographic, ethnic, or migratory factors, in
49 addition to pathogenesis and patient outcomes (8-11). Genotype A is highly prevalent across
50 Southeast Africa, Europe, North America, and India, whereas genotypes B and C are more commonly
51 found in the Asia-Pacific region. Genotype D is the most widespread and is prevalent in North Africa,
52 Europe, the Mediterranean region, North America, and India, while genotype E is primarily restricted
53 to West and Central Africa, as well as Saudi Arabia (12). Genotypes F and H are found in Central and
54 South America and are rare in other parts of the world (13), while the distribution of genotype G is
55 not completely known (14). Genotype I and J are found in Asia and probably result from recombination
56 events with other genotypes (15, 16).

57 **Table 1. Geographical distribution of HBV genotypes and subtypes.**58 *Adapted from Lin et.al, 2017 (10).*

Genotypes	Subtypes	Geographic location
A	A1 – A8	Sub-Saharan Africa, Asia, Northern Europe, America, India, Central and Western Africa, Gambia, Nigeria, Haiti
B	B1 – B10	East Asia, Taiwan, China, Indonesia, Vietnam, Philippines, Japan, Alaska, Northern Canada, Greenland
C	C1 – C17	Taiwan, China, Korea, Japan, Southeast Asia, Australia, Philippines, Vietnam, Indonesia, Papua New Guinea
D	D1 – D12	North Africa, Europe, Mediterranean countries, North America, India, Indonesia, Oceania
E		West and Central Africa, Saudi Arabia
F	F1 – F6	Central and South America
G		France, Germany, the United States
H		Central America
I	I1 – I3	Vietnam, Laos, southwestern China, eastern India
J		Japan

59

60 **HBV genotype F**

61 The S gene, preceded by a preS1 and S2 gene, produces the large, medium, and small forms of the
62 surface antigen (HBsAg). Genotype F displays marked differences from the other genotypes in pre-S/S
63 gene and studies have indicated that detecting genotype F may pose challenges due to the relatively
64 low sensitivity of assays designed for this particular type (1, 17). Genotype F is also the most divergent
65 of the genotypes and is subdivided into six subgenotypes (F1 to F6). Within subgenotypes F1 and F2,
66 distinct clades have been identified. As for subgenotype F1, multiple clades (F1a, F1b, F1c, and F1d)
67 have been described, each representing unique genetic lineages within the broader subgenotype
68 classification. Similarly, within subgenotype F2, distinct clades (F2a and F2b) have been identified.
69 These subclades likely reflect additional levels of genetic diversity and evolutionary divergence within
70 the overall framework of genotype F. (15). The phylogeny of genotype F is unique and shows little
71 intra-subgenotype diversity. However, long evolutionary distances are observed between its six

72 subgenotypes. This suggests that while the subgenotypes themselves exhibit relatively little genetic
73 variation, they have undergone significant evolutionary divergence from one another over time.

74 Moreover, there is evidence suggesting that genotype F might have the capacity to evade immunity
75 induced by vaccines (18). Vaccines against HBV primarily target the HBsAg and induce an immune
76 response that produces antibodies against HBsAg. However, one study has shown that HBV genotype
77 F strains may show variations in the HBsAg protein compared to other genotypes (19). These variations
78 could potentially affect the recognition and binding of antibodies produced by vaccinated individuals,
79 reducing the efficacy of vaccination against genotype F strains. However, more data is needed to
80 substantiate this hypothesis.

81 [HBV diagnosis](#)

82 The accurate detection of HBV infection is essential for early intervention, appropriate patient
83 management, and the prevention of transmission. In the Netherlands, HBV diagnosis is based on
84 detection of HBV antigen and antibodies in serum using serological immunoassays. Key markers, such
85 as HBV surface antigen (HBsAg), HBV envelope antigen (HBeAg), and antibodies to HBV core antigen
86 (anti-HBc), are predominantly tested to determine infection status.

87 The impact of HBV's high genetic variability on serological assay performance, specifically in detection
88 efficacy of HBsAg from different genotypes by assays used in routine diagnostics, remains an
89 understudied area.


90 [HBV surveillance in the Netherlands](#)

91 In the Netherlands, all acute HBV infections are reported anonymously by the GGD (Municipal or
92 Community Health Service) to OSIRIS (a public health database used for surveillance purposes).
93 However, clinical, patient, and test data on HBV markers are not reported, limiting the ability to
94 analyze the relationship between different HBV marker levels and substitutions. Acute cases are
95 reported to OSIRIS based on a positive HBsAg test result and/or an anti-HBc IgM result (if available).
96 Chronic cases are notifiable when a positive HBsAg or HBV DNA result is diagnosed for the first time
97 in the Netherlands. To determine risk exposure, trace source research and/or notify partners,
98 interviews are conducted with affected individuals. Since 2004, blood samples have been requested
99 for typing from all acute cases reported to OSIRIS. From 2010 onwards, samples from chronic cases
100 showing risk behavior (i.e., tested in the context of the HBV vaccination program for behavioral risk
101 groups or individuals reporting sexual contact by MSM as a route of transmission), were also
102 requested for typing (20).

103 In line with patterns observed in most European countries, genotype A is most prevalent in the
104 Netherlands, followed by genotype D (20, 21). Genotype A is consistently the most common genotype
105 each year, but the number of reported cases has generally decreased over the period from 77 cases
106 in 2011 to 50 cases in 2023 (Figure 1). Genotype D shows some fluctuation, but a noticeable decline
107 has also been observed in recent years, peaking at 30 cases in 2016, decreasing to 11 reported cases
108 in 2022, but saw an increase again to 25 cases in 2023. One of the most striking trends in the data is
109 the remarkable increase in the prevalence of genotype F, especially since 2019, even surpassing
110 genotype D as the second dominant type in 2021 (Figure 1). In the first part of the decade (2011-2018),
111 genotype F cases were relatively low and stable, typically ranging from 1 to 7 cases per year. As of
112 2019, there has been a significant increase in the number of reported cases of genotype F, with 13
113 cases in 2019, 7 in 2020, 14 in 2021, 9 in 2022, and 6 in 2023. This upward trend suggests that genotype
114 F is becoming increasingly common in the Netherlands.

115 This shift in genotype prevalence was identified in the Netherlands because of the implementation of
116 standard typing of all reported samples. Since 2004 blood samples from all acute cases are requested
117 for typing based on the S and C gene or complete genome (since 2017), facilitating the identification
118 of genotype variations and contributing to a more comprehensive understanding of the evolving
119 landscape of HBV infections in the region (Figure 1) (20).

120 **Figure 1. Annual distribution of HBV genotypes in the Netherlands over a period of 13 years,**
121 **between 2011 and 2023.**



122

123 Role of EQA/SKML

124 In the framework of promoting the quality of medical laboratory tests in the context of diagnosis and
125 treatment and to raise and maintain those tests at the highest possible level, the routinely designed
126 EQAs are distributed by the Foundation for Quality Assessment in Medical Laboratory Diagnostics
127 (Stichting Kwaliteitsbewaking Medische Laboratoriumdiagnostiek, SKML). EQAs of medical laboratory
128 diagnostics is required by ISO15189 as a means by which laboratories verify whether their methods
129 perform as expected and to what extent their results match those of other users of the same method.
130 The same data also provides insight into the performance of groups of users with the same method
131 and therefore into the performance of those methods. SKML is responsible for providing external
132 quality control distributions for all laboratory specialties, thereby promoting the standardization of
133 diagnostic procedures across laboratories. By participating in EQA programs, laboratories can align
134 their testing procedures with national and international standards, guaranteeing consistent and
135 reliable results.

136 Given the diversity and complexity of HBV genotypes, the need for a comprehensive HBV genotype
137 EQA panel has become increasingly evident. Traditional EQA panels may not adequately represent the
138 full genetic diversity of HBV, especially for less common genotypes. This can lead to gaps in diagnostic
139 accuracy and reliability, as assays may be less sensitive to these genotypes. The previous HBV
140 genotype panel used by SKML for Hepatitis B focuses on testing the measurement of HBsAg in serum,
141 with the positive/negative cut-off receiving a lot of attention, ensuring that laboratories can accurately
142 distinguish between positive and negative samples. However, the increase in the number of genotype
143 F cases in the Netherlands, has highlighted the possible need for more targeted assessment tools. By
144 incorporating a comprehensive HBV genotype panel, SKML can guarantee that laboratories can
145 accurately detect a wide range of HBV genotypes. This is crucial for effective patient management,
146 especially given the increasing diversity of HBV variants in the Netherlands.

147 Objectives

148 The surge in genotype F prevalence is reshaping the epidemiological landscape of HBV infections in
149 the Netherlands. The current diversification of genotypes in the Netherlands raises concerns about
150 potential vaccine evasion, especially with the influx of migrants from regions with a variety of
151 dominant genotypes.

152 In this study we aimed to investigate the serological detection performance on HBsAg diagnostic
153 assays used in Dutch Medical Microbiology and Clinical Chemical laboratories for different HBsAg
154 genotypes. For this we used the WHO HBV type specific panel (22) as a template for the design of an
155 EQA trial panel. The obtained EQA for different HBsAg genotypes will help in determining whether
156 diagnostic panels uniformly detect all HBV genotypes. In an ideal scenario, the assays should be able
157 to detect all included HBV genotypes equally well.

158

159 2. Study design and methods

160 Selection and characteristics of the EQA subgenotypes

161 An EQA panel consisting of sixteen well-characterized serological samples (15 HBV positive and 1
 162 negative serum control) with known HBV types and antigen amounts, representing various HBV
 163 subgenotypes was designed using WHO reference materials (PEI code 6100/09; version 3; Nov 13,
 164 2017) and published by Chudy et al. (1). The HBsAg concentration was determined by quantitative
 165 chemiluminescent immunoassay in IU/mL (CLIA, ARCHITECT HBsAg, Abbott, Germany). Detailed
 166 information on the panel is shown in Table 3. In detail, the materials used in the EQA consisted of non-
 167 inactivated freeze-dried plasma samples, obtained from the Paul-Ehrlich institute, a WHO
 168 collaborating centrum for quality assurance of blood products and *in vitro* diagnostic devices. As part
 169 of the blinding procedure, this information was not shared with the participating laboratories.

170 **Table 2. Characterization of the HBV genotype panel members and HBsAg measurements.**

171 *Target dilution to ca. 30 IU HBsAg/mL^a*

Sample	Origin	HBsAg (serological) subtype	HBV sub-genotype	HBsAg (IU/mL) ^a
1	Brazil	adw2	A1	29.23
2	Germany	ayw2	D1	27.77
3	Germany	adw2	A2	26.77
4	West Africa	ayw4	E	28.50
5	Japan	adw2	B2	24.90
6	South Africa	adw2	A1	31.80
7	Japan	adr	C2	28.40
8	Russia	adr	C2	30.00
9	South Africa	ayw2	D3	24.13
10	None	/	/	0.00 (NEG CTRL)
11	Japan	adr	C2	29.23
12	Brazil	adw4	F2	31.67
13	Russia	ayw3	D2	29.40

14	Japan	adw2	B2	21.10
15	Brazil	adw4	F2	26.80
16	Germany	adw4	H	29.50

172 ^a1 in 3 dilutions based on mean CLIA values in IU/mL (ARCHITECT QT, Abbott)(23)

173

174 Procedure for participating laboratories

175 Identification of eligible laboratories

176 This EQA was set up as a pilot project in collaboration with SKML and Medisch Centrum Alkmaar
177 (MCA). We distributed invitations with forms for collection data on available testing systems;
178 Appendix 1) to laboratories to perform serological testing on samples containing different HBV
179 genotypes. All laboratories expressing interest were included to consider diverse detections systems
180 (both quantitative and qualitative assays). Both the participating laboratories as the reported assays
181 received a specific laboratory and test code.

182 Distribution of EQA subgenotypes to laboratories and reporting of results

183 All samples were first allowed to equilibrate to ambient temperature, before being resuspended in
184 0.4 mL of sterile nuclease-free water. This dilution step ensured that the assays were challenged near
185 their lower detection limits, providing a more rigorous evaluation of their performance. Samples were
186 then left for a minimum of 10 minutes with occasional agitation before being aliquoted into new,
187 relabeled vials (maintaining a new, randomized sample order). All samples were stored and shipped
188 as you would with handling patient samples.

189 Upon arrival, diagnostic laboratories were requested to check all samples for damage, completeness,
190 and shipping conditions, and to store the samples like patient material until testing. The samples were
191 deemed stable up to 5 days after reception if stored under defined conditions (i.e., refrigerated).
192 Laboratories were advised to freeze samples when processing was expected to start only after five
193 days. MCA was responsible for the distribution of the EQA panel to the participating laboratories.

194 The distributed panel was accompanied by detailed instructions for testing, consisting of data
195 collection forms (Appendix 2) and a testing protocol (Appendix 3). The diagnostic laboratories
196 participating in the EQA were instructed to assess the panels using each of their routine serological
197 procedures and to report the results obtained for each of the sixteen samples using the provided form
198 (Appendix 3). Results should be reported in IU/mL for quantitative tests and in sample-to-cutoff (S/CO)
199 ratios for qualitative tests. Data about the type of assay utilized, detection methods, instrument and

200 reagent details, manufacturers, raw detection data, quantitative, and qualitative results of HBsAg
201 were requested as well.

202 Analysis of the results

203 The data collection was performed through the RIVM during 2023 up to January 2024 and the data
204 obtained from each diagnostic laboratory were reported as total number of participating laboratories
205 reporting data. The amount of different serological antigenic assays was also reported (Table 3). Per
206 each serological assay the results were analyzed.

207 Fictional standardization

208 The results obtained were further analyzed through intra-assay fictional standardization to identify
209 the causes of differences in results. The idea behind fictional standardization is that if each laboratory
210 had included the same calibrator-sample in the measurements, the results could have been converted
211 to that calibrator-sample. If those converted values were then remarkably close to each other, the
212 original difference between the results of the participating laboratories could be explained by
213 difference in standardization. This can also be done fictionally with the obtained results. Let say that
214 we have n participating laboratories and each laboratory measures a value m_{ij} (the individual
215 measurement from laboratory i for sample j). If we had made a mixture with equal parts of all tested
216 samples, and measured this mixture by all participating laboratories, the mixture value (MV) found
217 would have been equal to the average of all the measured values. This essentially normalizes the
218 measurements across different laboratories, allowing for a more direct comparison. It is like creating
219 a “composite sample” that represents the average of all samples. In this way, the measurement from
220 each individual laboratory can be directly compared to this “composite sample” measurement.

221 In this study, the average value of all measured subgenotype samples across all laboratories (a_g) was
222 calculated. From there, the average value of all measured values across all subgenotypes was
223 calculated by dividing a_g by the number of samples (s) in our panel (i.e., 16 – 1 negative control = 15).

$$224 \quad V = \frac{a_g}{s}$$

225 This value (V) is then the assigned value of the fictitious calibrator.

226 Now, to calculate the calibration factor c_i for each laboratory, the fictitious calibrator (V) is divided by
227 a_i , where a_i is the average value measured per each laboratory (i) across all subgenotypes.

$$228 \quad c_i = \frac{V}{a_i}$$

229 This calibration factor c_i allows us to standardize the values reported by different laboratories.

230 As explained above, m_{ij} is the individual measurement from laboratory i for sample j . This
231 measurement is then multiplied by the calculated calibration factor c_i to get the fictitious calibrated
232 value f_{ij} for sample j from laboratory i .

233
$$f_{ij} = c_i \cdot m_{ij}$$

234 In this way, all measurements are normalized and can be directly compared across the different
235 participating laboratories.

236

237


3. Results

238 In 2023, the EQA had been tested by 27 laboratories across the Netherlands, including the Caribbean
239 islands Aruba, Curaçao, and Sint-Maarten (Figure 2). The results were evaluated for all these
240 laboratories. The HBV subgenotype panel for HBsAg contains 16 samples representing the most
241 prevalent HBV subgenotypes (Table 1). The HBV DNA concentration in the samples ranged from 21.10
242 (sample 14) to 31.80 IU/mL (sample 6), with a median of 28.50 IU/mL.

243 **Figure 2. Numbering of participating MMLs returning external quality assessment results, by**
244 **municipality.**

Participating laboratories in the Netherlands

Aruba, Sint-Maarten, and Curaçao not shown



245

246 Table 3 and 4 report the number of all the different diagnostic test assays used by the participating
247 laboratories and the units of reported results.

248 A wide range of assays was used across the participating laboratories (Table 3), with laboratories
249 reporting results in different units, such as S/CO ratio, IU/mL, and reflective lights unit (RLU) (Table 4).
250 In total, 27 laboratories submitted 13 quantitative data sets obtained with 3 different quantitative
251 tests and 21 qualitative data sets from 8 different qualitative tests. One laboratory (HBVg18)
252 performed a quantitative HBsAg test (Liaison-XL Murex HBsAg Quant, test code 2A and 2B) and
253 reported results in both IU/mL and in reflective light units (RLU). Laboratories HBVg06, 11, 14, 17, 19,
254 22, 23, and 25 reported results in multiple units for the same assay, such as both RLU and S/CO ratio.
255 Four laboratories (HBVg02, 09, 17, and 22) used multiple assays and reported in several units. Dataset
256 from 5 laboratories (HBVg09, 15, 25, 11, and 14) were excluded from analysis as they did not provide
257 S/CO ratio values when a single measurement exceeded 1000 S/CO ratio. Despite this, S/CO ratios \geq
258 1000 were considered positive. These datasets were reported from using either the Advia Centaur XPT
259 HBsAg II (test code 6B, 6C) or the Atellica IM analyzer HBsAg II (test code 6A, 6D) assays (both from
260 Siemens Healthcare Diagnostics inc.).

261

262 **Table 3. List of all assay kits used by the participating MMLs.**

Test no	HBsAg test	Kit manufacturer	Test code
Quantitative assays			
1	Liaison-XL Murex HBsAg Quant	DiaSorin S.p.A.	2A
2	Elecsys HBsAg II Quant using Cobas 8000 modular analyzer systems	Roche Diagnostics Division	3A
3	Alinity-i HBsAg II QT	Abbott Diagnostics Division	4A
Qualitative			
4	Advia Centaur XPT HBsAg II	Siemens Healthcare Diagnostics Inc.	6B, 6C
5	Atellica IM Analyzer HBsAg II	Siemens Healthcare Diagnostics Inc.	6A, 6D
6	Vidas HBsAg ULTRA	bioMérieux	5
7	ARCHITECT HBsAg II Qual	Abbott Diagnostics Division	4B
8	Alinity-i HBsAg II Qual	Abbott Diagnostics Division	4C
9	Liaison-XL Murex HBsAg Qual	DiaSorin S.p.A.	2B
10	UniCel DxI 600 Access 2 HBsAg	Beckman Coulter Diagnostics	1
11	Elecsys HBsAg II Qual using Cobas 6000 or 8000 modular analyzer systems	Roche Diagnostics Division	3B, 3C

263

264

265

266 **Table 4. Reported assays used for HBV genotype testing by all participating laboratories.**

Lab code	Used assay (by test code)	Results reported in
HBVg01	1	S/CO ratio
HBVg02	2A	IU/mL
	1	S/CO ratio
HBVg03	2A	IU/mL
HBVg04	3C	S/CO ratio
HBVg05	2A	IU/mL
HBVg06	5	RLU ^a
		S/CO ratio
HBVg07	4B	S/CO ratio
HBVg08	2A	IU/mL
	3C	S/CO ratio
HBVg09	6A/D	S/CO ratio
	2A	IU/mL
HBVg10	2A	IU/mL
		RLU ^a
HBVg11	6B/C	S/CO ratio
		IU/mL
HBVg12	2A	IU/mL
		IU/mL
HBVg13	2A	IU/mL
		RLU ^a
HBVg14	6B/C	S/CO ratio
		S/CO ratio
HBVg15	6A/D	S/CO ratio
HBVg16	3C	S/CO ratio
HBVg17	4A	IU/mL
	4C	RLU ^a

		S/CO ratio
HBVg18	2A	IU/mL
	2B	RLU ^a
HBVg19	4C	RLU ^a S/CO ratio
	3C	S/CO ratio
HBVg20	3C	S/CO ratio
HBVg21	2A	IU/mL
HBVg22	3A	IU/mL
	3C	S/CO ratio
	4B	RLU ^a S/CO ratio
HBVg23	4C	RLU ^a S/CO ratio
	5	RLU ^a S/CO ratio
HBVg24	2A	IU/mL
HBVg25	6A/D	RLU ^a S/CO ratio
HBVg26	3B	S/CO ratio
HBVg27	2A	IU/mL

267 ^aRLU: Reflective Lights Unit

268





269 Fifteen qualitative data sets were analyzed (Supplementary Figure 1), with the Elecsys HBsAg II Qual
 270 assay (test code 3B, 3C; Roche Diagnostics) being the most frequently used assay and reported by six
 271 different laboratories (Table 4). As depicted in Figure 3A, consistent detection of all HBV genotypes,
 272 including genotype F2, was demonstrated by all assays. Of interest, the Vidas HBsAg ULTRA assay (test
 273 code 5) reported S/CO ratios at a much lower magnitude (average S/CO ratio of 24.70) compared to
 274 the other four assays. Excluding the Vidas HBsAg ULTRA assay (test code 5), the lowest S/CO ratio was

275 obtained when using the UniCel Dxl 600 Access assay (test code 1) and was lowest for sample 2
276 (127.72; HBV subgenotype D1). Additionally, both ARCHITECT HBsAg II (test code 4B) and Alinity-I
277 HBsAg (test code 4C; Abbott) yielded the highest S/CO ratios for the tested samples. Notably, the
278 Alinity-I HBsAg II QL assay (test code 4C) exhibited the highest standard deviation (SD) in reported
279 values.

280 Three quantitative assays, reporting 13 data sets, were analyzed (Supplementary Figure 2). The
281 Liaison-XL MUREX assay (test code 2A) was by far the most frequently reported assay, with as many
282 as eleven laboratories utilizing this specific assay (Table 4). Conversely, the other two assays were each
283 reported only once. Consistent detection of all HBV genotypes was observed across all assays (Figure
284 3B). The Liaison-XL Murex quantitative assay (test code 2A) displayed variability in antigen detection
285 for the same subgenotype among laboratories employing this specific diagnostic system.
286 Furthermore, all quantitative assays exhibited variations in HBsAg detection rates, with values ranging
287 from 10 (subgenotype A2) to 76 (subgenotype E) IU/mL. When using the Liaison-XL Murex quantitative
288 assay (test code 2A), the lowest measured values (average HBsAg < 29 IU/mL) were observed for
289 subgenotypes A2, B2, D1, D2, D3, and F2, whereas the highest values were noted for subgenotypes
290 C2 with Japanese origin (average HBsAg of 51 IU/mL) and E with West African origin (average HBsAg
291 of 45 IU/mL).

292 To distinguish whether the variation between genotypes originated from variations in intra-assay
293 standardization procedures among participating laboratories or from differences in the type of sample
294 measured, we applied the concept of fictional standardization (see methods section for detailed
295 explanation). Fictional standardization was applied to both qualitative and quantitative assays, and as
296 demonstrated in Figure 3 C and D, after fictional standardization, the SD was reduced compared to
297 data presented in Figure 3 A and B, thus reducing the variability across the different assays used. For
298 instance, HBV subgenotype A1 had an average IU/mL value of 26.91 across all laboratories using
299 Liaison-XL Murex (test code 2A), with an initial SD of 5.9. Following fictional standardization, the SD
300 decreased to 3.7.

Figure 3. Comparative analysis of qualitative (A and C) and quantitative (B and D) diagnostic HBsAg serological assays, before (A and B) and after (C and D) fictional standardization. For A and C, HBsAg concentration is expressed as S/CO ratio. For B and D, HBsAg quantity is shown in IU/mL. Mean values \pm standard deviation (SD) is shown. In brackets the number (n) indicates the total number of datasets reported per assay. In B and D, the reference HBsAg values (Table 2), adjusted from Chudy et al. (1) provide a benchmark for comparison. Note: Test code 5 (Vidas HBsAg ULTRA assay) displays significantly lower S/CO ratios, appearing near zero despite actual non-zero values.



Observed challenges

Limitations in reporting S/CO ratio values

One of the challenges observed was the limitations in reporting S/CO ratio values. The sample-to-cutoff ratio is a measure used in immunoassays to determine the presence of an analyte. This is calculated by dividing the signal of the test sample by the signal of a defined internal control (cutoff). A ratio above a certain threshold indicates a positive test result.

In our study, five laboratories reported the use of either the Advia Centaur XPT HBsAg II or the Atellica IM analyzer HBsAg II assays (both Siemens Healthcare Diagnostics inc., respective teste code number 6B/6C and 6A/6D). These five data sets were excluded as they did not specify S/CO ratio values when a single measurement exceeded 1000 S/CO ratio. This limitation hindered our ability to evaluate and compare the performance of these assays with others, particularly in instances of high viral load.

Inconsistent assay usage

An inconsistency was noted in the usage of the Liaison-XL Murex HBsAg by one laboratory for qualitative reporting. This test is primarily used as a quantitative assay. Because of this deviation from standard practice, the results from this laboratory using this assay for qualitative reporting were excluded in the analysis.

Communication obstacles

One laboratory became unreachable for further inquiries or comments regarding their submitted report. This communication barrier hindered the clarification of certain aspects of their data.

4. Conclusions

EQA programs play a crucial role in evaluating the performance and status of diagnostic assays in clinical laboratories (24). Ensuring the reproducibility and reliability of diagnostic assays are of particular importance for effective clinical management and public health initiatives (24). Laboratories should conduct verification studies before and continuously monitor reliability through routinely quality management processes (25). Moreover, the inter-laboratory comparison facilitates the assessment of participants' performance, enabling the evaluation of analytical performance, test interpretation, and method performance. By implementing these measures, laboratories can enhance the accuracy and consistency of diagnostic testing.

The results of the EQA presented here was organized in 2023 to evaluate the serological detection performance on HBsAg assays used for different genotypes in routine diagnostics in Dutch Medical Microbiology and Clinical Chemical laboratories. This HBV genotype panel was specifically designed to represent the wild type strains of the most prevalent HBV subgenotypes. The 27 laboratories participating in this EQA program were required to test the panel using their routine procedures and report their qualitative and/or quantitative results. A total of eleven different diagnostic assays were evaluated, three quantitative and eight qualitative assays, covering 34 data sets. The EQA demonstrated consistent and adequate serological detection of Hepatitis B surface antigen from all tested HBV genotypes, including F2, across the wide range of assays used. The results from this study indicate good assay performances of HBsAg immunoassays. However, as demonstrated in Figure 3 C and D, after fictional standardization, the SD was reduced compared to data presented in left panels of the same figure (Figure 3 A and B). This indicates that the discrepancies observed in the reporting results can be attributed, at least partly, to the utilization of various standardization methods or conventional instrumental-specific calibrators across participating laboratories. Variations in the calibration protocols, reference materials, or calibration intervals used by different laboratories may introduce differences in assay performance. In addition, differences in instrument models and/or manufacturers may require adjustments in calibration procedures, further contributing to these discrepancies in the observed results. It is noteworthy that a larger SD is observed when considering variations between different genotypes. This is more pronounced for subgenotypes such as C2, or subgenotype E, where greater intra- and inter-variability is observed between different test systems compared to other subgenotypes.

5. Recommendations

Addressing the challenges observed is essential for improving the overall quality and consistency of diagnostic practices for HBV in Dutch medical microbiology laboratories. Several strategies can be implemented to achieve this objective. Standardization of reporting is one option. For example, working closely with the assay manufacturers to establish clear guidelines for reporting S/CO ratio values, particularly for measurements exceeding 1000. This may include updating assay protocols or developing additional guidelines to ensure consistent reporting practices across laboratories.

Another option is to invest in capacity-building with harmonized practices. Provide additional training and educational sessions for laboratory personnel to ensure proficiency in the appropriate use and interpretation of diagnostic assays. This would lead to minimizing discrepancies in the reporting of assays.

Having clear communication protocols between medical microbiology and clinical chemical laboratories on the one hand and EQA organizers on the other, would help ensure timely and effective information exchange. One option is to use an online webtool or platform for submitting results, supplemented with an electronic feedback survey for participating laboratories, instead of using a preformatted Excel sheet.

Another recommendation is the development of follow-up procedures for laboratories that become unreachable or fail to provide adequate information during the EQA process. Providing assistance to laboratories experiencing difficulties is one option. This could include technical support, guidance on what to do, on how to submit required information, or addressing other obstacles they may face. This would enable proactive resolution of encountered barriers. By implementing these strategies, the EQA program can improve the proficiency and consistency of diagnostic practices among participating laboratories.

Consistent HBV diagnostics is especially important when considering antiviral testing. While HBsAg serology levels provide valuable information, the decision to initiate antiviral therapy or assess the efficacy of antiviral treatments primarily depends on accurate measurements of HBV DNA levels (26, 27). HBV DNA quantification is the key marker used to assess viral replication and determine the need for antiviral therapy. Discrepancies in assay results underscore the need for improved standardization or calibration methods across different laboratories and assay systems. Accurate calibration ensures that all laboratories report consistent and comparable results, which is crucial for making informed treatment decisions. Using data from these and similar EQAs to inform and enhance HBV diagnostics

will improve the reliability of antiviral testing schemes and assist in making reliable decisions about initiating antiviral therapy.

Acknowledgments

We thank Dr. C. Weykamp for help with the analysis of the fictional standardization. Dr. G. Boland for being the responsible partner and contact person for HBV serological EQA for SKML. We thank all partners and technicians of the participating medical microbiology laboratories involved in this study for the consent to public the results of this EQA study. A special acknowledgment to Dr. H. ten Bruin from MCA for all the help and the coordination of the study.

Disclaimer

The author is a fellow of the ECDC Fellowship Programme, supported financially by the European Centre for Disease Prevention and Control (ECDC). The views and opinions expressed herein do not state or reflect those of ECDC. ECDC is not responsible for the data and information collation and analysis and cannot be held liable for conclusions or opinions drawn.

References

1. Chudy M, Scheiblauer H, Hanschmann KM, Kress J, Nick S, Wend U, et al. Performance of hepatitis B surface antigen tests with the first WHO international hepatitis B virus genotype reference panel. *J Clin Virol*. 2013;58(1):47-53.
2. Organization WH. Hepatitis B, factsheet 2024 [Available from: <https://www.who.int/news-room/fact-sheets/detail/hepatitis-b>].
3. Soriano V, Moreno-Torres V, Trevino A, de Jesus F, Corral O, de Mendoza C. Prospects for Controlling Hepatitis B Globally. *Pathogens*. 2024;13(4).
4. Jeng WJ, Papatheodoridis GV, Lok ASF. Hepatitis B. *Lancet*. 2023;401(10381):1039-52.
5. Zhang ZH, Wu CC, Chen XW, Li X, Li J, Lu MJ. Genetic variation of hepatitis B virus and its significance for pathogenesis. *World J Gastroenterol*. 2016;22(1):126-44.
6. Kramvis A. Genotypes and genetic variability of hepatitis B virus. *Intervirol*. 2014;57(3-4):141-50.
7. Gerlich WH, Glebe D, Kramvis A, Magnius LO. Peculiarities in the designations of hepatitis B virus genes, their products, and their antigenic specificities: a potential source of misunderstandings. *Virus Genes*. 2020;56(2):109-19.
8. Lin CL, Kao JH. Hepatitis B virus genotypes and variants. *Cold Spring Harb Perspect Med*. 2015;5(5):a021436.
9. Kao JH. Hepatitis B viral genotypes: clinical relevance and molecular characteristics. *J Gastroenterol Hepatol*. 2002;17(6):643-50.
10. Lin CL, Kao JH. Natural history of acute and chronic hepatitis B: The role of HBV genotypes and mutants. *Best Pract Res Clin Gastroenterol*. 2017;31(3):249-55.
11. Liu Z, Zhang Y, Xu M, Li X, Zhang Z. Distribution of hepatitis B virus genotypes and subgenotypes: A meta-analysis. *Medicine*. 2021;100(50):e27941.
12. Ingasia LAO, Wose Kinge C, Kramvis A. Genotype E: The neglected genotype of hepatitis B virus. *World J Hepatol*. 2021;13(12):1875-91.
13. Mojsiejczuk LN, Torres C, Pisano MB, Re V, Campos RH, Flichman DM. New pieces on genetic diversity and evolutionary history of hepatitis B virus: Characterization of the novel subgenotype F6. *Infect Genet Evol*. 2017;47:140-2.
14. Wolf JM, De Carli S, Pereira V, Simon D, Lunge VR. Temporal evolution and global spread of hepatitis B virus genotype G. *J Viral Hepat*. 2021;28(2):393-9.
15. Toye RM, Loureiro CL, Jaspe RC, Zoulim F, Pujol FH, Chemin I. The Hepatitis B Virus Genotypes E to J: The Overlooked Genotypes. *Microorganisms*. 2023;11(8).
16. Tatematsu K, Tanaka Y, Kurbanov F, Sugauchi F, Mano S, Maeshiro T, et al. A genetic variant of hepatitis B virus divergent from known human and ape genotypes isolated from a Japanese patient and provisionally assigned to new genotype J. *J Virol*. 2009;83(20):10538-47.
17. Limeres MJ, Gomez ER, Nosedá DG, Cerrudo CS, Ghiringhelli PD, Nusblat AD, Cuestas ML. Impact of hepatitis B virus genotype F on in vitro diagnosis: detection efficiency of HBsAg from Amerindian subgenotypes F1b and F4. *Arch Virol*. 2019;164(9):2297-307.

18. Schlabe S, van Bremen K, Aldabbagh S, Glebe D, Bremer CM, Marsen T, et al. Hepatitis B virus subgenotype F3 reactivation with vaccine escape mutations: A case report and review of the literature. *World J Hepatol.* 2018;10(7):509-16.
19. Di Lello FA, Ridruejo E, Martinez AP, Perez PS, Campos RH, Flichman DM. Molecular epidemiology of hepatitis B virus mutants associated with vaccine escape, drug resistance and diagnosis failure. *J Viral Hepat.* 2019;26(5):552-60.
20. Cremer J, van Heiningen F, Veldhuijzen I, Benschop K. Characterization of Hepatitis B virus based complete genome analysis improves molecular surveillance and enables identification of a recombinant C/D strain in the Netherlands. *Heliyon.* 2023;9(11):e22358.
21. Hahné S, van Houdt R, Koedijk F, van Ballegooijen M, Cremer J, Bruisten S, et al. Selective Hepatitis B Virus Vaccination Has Reduced Hepatitis B Virus Transmission in The Netherlands. *PLOS ONE.* 2013;8(7):e67866.
22. Paul-Ehrlich-Institut. First WHO international reference panel for HBV genotypes for Hepatitis B surface Antigen assays. 2017.
23. Chudy M, Hanschmann KM, Kress J, Nick S, Campos R, Wend U, et al. First WHO International Reference Panel containing hepatitis B virus genotypes A-G for assays of the viral DNA. *J Clin Virol.* 2012;55(4):303-9.
24. Ceriotti F. The role of External Quality Assessment Schemes in monitoring and improving the standardization process. *Clin Chim Acta.* 2014;432:77-81.
25. Lisboa Bastos M, Tavaziva G, Abidi SK, Campbell JR, Haraoui LP, Johnston JC, et al. Diagnostic accuracy of serological tests for covid-19: systematic review and meta-analysis. *BMJ.* 2020;370:m2516.
26. WHO. GUIDELINES FOR THE PREVENTION, CARE AND TREATMENT OF PERSONS WITH CHRONIC HEPATITIS B INFECTION. 2015.
27. Guvenir M, Arikan A. Hepatitis B Virus: From Diagnosis to Treatment. *Pol J Microbiol.* 2020;69(4):391-9.

Appendices

Appendix 1. Invitation letter.

Call for interest – Serologic hepatitis B EQA: assessment of diagnostic capacity for distinct genotypes

Dear colleagues,

Since 2019, an uprise of HBV genotype F has been noted in the Netherlands (and no data from other countries is available to our knowledge). Most recently, this genotype detection surpassed the most common genotype D, giving place to a new epidemiological picture of HBV infections which we would like to better understand.

The aim is to study the detection capacity of distinct HBV genotypes and subgenotypes by using different diagnostic systems. At the same time, this EQA will allow the collection of valuable information of ongoing monitoring/surveillance activities and foster preparedness (capacity to detect uncommon genotypes).

We would like to ask for your interest to collaborate in testing your HBsAg detection assay with an HBV serological panel.

Unfortunately, as funds are limited, we will only be able to distribute the panel to 12-15 laboratories expressing interest. In case more than 15 laboratories have expressed interest we will have to make a selection, taking into account diverse detection systems (quantitative assays, preferential). The panel (maximum of 18 samples) will be sent out by the end of July with a testing period of 3 weeks.

We would appreciate if you could complete the attached short survey to indicate whether you would like to participate or not. Please return the survey to us **before the 23rd June 2023**, to kim.benschop@rivm.nl and maggie.pires.simoës@rivm.nl.

A final decision will be shared among those willing to contribute by end of June.

This study is performed in collaboration with SKML. Results will be collected, analysed, and reported via RIVM.

Thank you in advance for your collaboration.

Yours sincerely,

Kim Benschop, Maggie Simões

Laboratory details

Name of the laboratory: _____

City / Province: _____

Contact name: _____

E-mail, phone number: _____

Interested in participating? **Please check your choice:** Yes ____, No ____*If not interest to participate in this Ring Trial, no further answers are needed.*

Diagnostic method	Diagnostic reagents/assay	System (Manufacturer and Model)	Min. sample volume input required (in μL)	Detection limits (IU; ng/mL)	Quantitative assay?	Lowest reliable detection limit	Assay cut-off to define positive and negative result	HBV Ag sample matrix validated (manufacturer, own validation)
					Yes / No			

E-mail this form to: kim.benshop@rivm.nl & maggie.pires.simoese@rivm.nl**Thank you for your collaboration**

Appendix 2. Reporting form.

Results of Ring Trial / Laboratory Proficiency Testing

Hepatitis B subgenotype serologic detection

Laboratory code	
Lab name	
Responsible name	
Responsible email	
Diagnostic System/Method	
Diagnostic reagents/ assay (kit)	
Diagnostic system (Manufacturer and Model)	
Type of assay please select from drop down list	If Qualitative assessment, further reporting in row 17 If Quantitative assessment, further reporting in row 40
Panel received (date: dd-mm-yyyy)	
Panel condition	If issues, please specify here:
Testing period (date: dd-mm-yyyy to dd-mm-yyyy)	

Qualitative with cutt-off							
Sample number	OD sample	OD cut-off	S/CO ratio	Results		Result/Interpretation	Remark if applicable. Clarify if technical issues
				HBsAg concentration (value)	Units		
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							

NOTES (if needed)	
Disclaim any differences to the manufacturer's instructions	
General observations	

NB: Columns in yellow consisted of drop-down choices to be selected from.

Quantitative					
Sample number	Results		Result / Interpretation	Remark if applicable	Clarify if technical issues
	HBsAg concentration (value)	Units			
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					

NOTES (if needed)	
Disclaim any differences to the manufacturer's instructions	
General observations	

NB: Columns in yellow consisted of drop-down choices to be selected from.

Appendix 3. Instructions for use of the trial EQA panel.

Please find enclosed the samples for the Pilot serologic HBV EQA panel

Contact Information:

For all questions please contact: Kim Benschop and Maggie Simoes
kim.benschop@rivm.nl; maggie.pires.simoaes@rivm.nl; +31 629637711

Control:

Please check all samples for damage, completeness, and shipping conditions.

These samples were sent without specified shipping conditions.

Special Comments:

An email will be sent with your assigned laboratory code and reporting template.
All communications should have your unique identifier (Lab Code) in the subject.

Closing date results:

Results are to be reported until the 15th.September.2023 (end of working day).

Responsible coordinator:

Dr. Kimberley S.M. Benschop
Centre for infectious Disease Control (Cib) | Centre for Infectious Disease Research, Diagnostics, and laboratory Surveillance (IDS).
Dutch National Institute for Public Health and the Environment
PO Box 1 | 3720 BA Bilthoven
Antonie van Leeuwenhoeklaan 9 | 3721 MA Bilthoven

Contents of this shipment:

One box labeled "serologic HBV EQA panel, Samples 1 to 8" and one box "serologic HBV EQA panel, Samples 9 to 16".

Each box contains 8 vials with each 0.4 mL of sample material.

Storage:

The samples must be stored like patient material.

The shelf life of the samples equals patient material.

Expiry date samples:

The samples are stable up to 5 days of reception if stored under defined conditions.

Preparation for use:

- a. Mix carefully.
- b. Process as patient material.

Measurement of samples:

The samples should be treated in the same manner as regular patient samples and in accordance with the regular used method(s).

Submission of results:

Results should be sent to : kim.benschop@rivm.nl; maggie.pires.simoaes@rivm.nl until the 15th.September.2023 (end of working day)

Risks:

Like patient specimens the material should be regarded as potentially infectious and be treated as such.

Appendix 4. Steps followed for fictional standardization

1. **Collect Data:**

Gather measurement results from each laboratory for the samples of interest.

2. **Calculate Averages:**

Compute the average value for each sample within each laboratory, and then calculate the average of these values across all labs. This gives you the assigned value of the fictitious calibrator.

3. **Apply Calibration Factor:**

Divide each laboratory's average value by the assigned value of the fictitious calibrator. This gives you a calibration factor for each laboratory.

4. **Apply Calibration Factor to Individual Measurements:**

Multiply each individual measurement from each laboratory by the corresponding calibration factor obtained in step 3.

5. **Assess Variability:**


Calculate the standard deviation (SD) of the measurements before and after calibration. Compare the SDs to assess the impact of standardization differences on the variability of the measurements.

6. **Analyze Results:**

Interpret the changes in variability and assess the extent to which differences in standardization contribute to the variability in the measurements.

Supplementary Figure 1.

Individual graphs of the 5 qualitative HBsAg assays used to test the 15 serological HBV genotype samples (excluding the negative control). HBsAg concentration is expressed as S/CO ratio.



Supplementary Figure 2.

Individual graphs of the 3 quantitative HBsAg assays used to test the 15 serological HBV genotype samples (excluding the negative control). HBsAg quantity is shown in IU/mL.

