

Moleculaire diagnostiek van helminthen: een pilot rondzending voor detectie in feces

20 Maart 2018

Lisette van Lieshout

namens de sectie parasitologie van de SKML



NEGLECTED TROPICAL DISEASES
SUPPORT CENTER



Stop Anthelmintic
Resistant Worms

Piet Cools &
Bruno Levecke



**GHENT
UNIVERSITY**

Jaco Verweij

Bestaande EQAS voor feces protozoa PCR

Clin Chem Lab Med 2018; aop

Theo A. Schuurs ^a, Rob Koelewijn, Eric A.T. Brienen ^a, Titia Kortbeek ^a, Theo G. Mank ^a, Bert Mulder ^a, Foekje F. Stelma ^a, Lisette van Lieshout ^a and Jaap J. van Hellemond ^{a,*}

Harmonization of PCR-based detection of intestinal pathogens: experiences from the Dutch external quality assessment scheme on molecular diagnosis of protozoa in stool samples




Figure 3: Reported Cq-values for *Giardia lamblia* in stool (black dots) and purified DNA (blue dots) specimens.

The reported Cq-values are plotted on the y-axis expressed as the reported Cq-value of the individual participant subtracted by the median of reported Cq-value of all participants. Each lane represents a single participant and the order in which the participants are plotted on the x-axis differs from that in Figures 1 and 2.

Waarom een EQAS voor wormen PCR?

Voordelen microscopisch feces onderzoek wormen:

- Brede diagnostiek, zoveel verschillende species

Nadelen microscopisch onderzoek wormen:

- Komt relatief erg weinig voor in NL
- Hoe behouden van competentie
- €€€€



Waarom een EQAS voor wormen PCR?

Nadelen PCR

- Teveel potentiële targets (nooit compleet)
- €€€€; zeker indien nooit positief
- (Nog) Niet 24/7 of zelfs 8/5 beschikbaar
- Beperkt aantal centra heeft momenteel de capaciteit/kundigheid

Welke targets hoogste prioriteit?

- Welke in gebruik?
- Welke behoefte aan EQAS?



Welke helminth targets in een feces PCR?

Klinische setting?

- Welke het makkelijkst gemist?
- Welke klinisch het meest relevant?

Survey setting?

- “Identification of hot-spots” & “post-MDA monitoring”
- Welke het makkelijkst gemist?
- Welke het meest relevant?

Internationale markt voor een EQAS?

- Specifieke logistieke uitdagingen?
- Klinische materialen?

Is er behoefte aan helminthen PCR Klinische setting in NL? Welk target eerst?

Transactions of the Royal Society of Tropical Medicine and Hygiene (2009) 103, 967–972



ELSEVIER

available at www.sciencedirect.com



journal homepage: <http://www.elsevier.com/locate/trstmh>



REVIEW

Strongyloidiasis – the most neglected of the neglected tropical diseases?


Annette Olsen^{a,*}, Lisette van Lieshout^b, Hanspeter Marti^c, Ton Polderman^b,
Katja Polman^d, Peter Steinmann^{e,f}, Russell Stothard^g, Søren Thybo^h,
Jaco J. Verweij^b, Pascal Magnussen^a

Diagnosing *Strongyloides* in a Dutch setting

Serology

- Different formats (in-house, commercial)
- Screening of specific patients, migrants, chronic infections

Antibody detection



Microscopy?


- Travelers (based on LUMC serology)
- Post-treatment monitoring (disseminated infections)

Diagnosis of *Strongyloides*

Microscopy based detection of larvae (no eggs!)



sensitivity stool culture,
N=19 patients



Nielsen & Mojon 1987

Differentiation between *Strongyloides* L1 and L3 larvae

L1 = rhabditoid larvae



lengths: 200-300 μm
short buccal cavity
large genital primordium

L3 = filariform larvae




lengths: 500-600 μm
motile, slender
long oesophagus ($>1/3$)
no sheath, notched tail

Peru: *Strongyloides* survey

3 consecutive stool samples

PCR: Verweij *et al* (2009) – 18S ribosomal RNA gen



* Baermann

La Merced
(N=188)


Verweij *et al* (un published data)

Mozambique: *Strongyloides* survey

PCR > microscopy, but not 100% sensitivity

Meurs et al., 2017 Plos NTD

Strongyloides PCR: Verweij *et al* (2009) – 18S ribosomal RNA gen



Results Antwerp Travel Clinic

N=2591 complete data

	Microscopy	PCR
<i>E. histolytica/E. dispar</i>	99	
<i>E. histolytica</i>		13
<i>Giardia lamblia</i>	95	149
<i>Cryptosporidium</i>	12	31
<i>Strongyloides stercoralis</i>	3*	21**

*) Baermann performed in 121 clinically suspected cases only

***) positive in microscopy N=3, serology N=7. Cases with eosinophilia N=5, clinical presentation N=4

Diagnosing *Schistosoma* in a Dutch setting

Serology

- Different formats (in-house, commercial)
- Screening of specific patients, travelers


Antibody detection



Microscopy?


- Migrants, chronic infection (based on LUMC serology)
- Post-treatment monitoring

S. haematobium school children Ghana (N=730)



Microscopy: eggs detected in 10 mL of urine


S. haematobium school children Ghana (N=730)



Schistosoma DNA isolated from 200 μ l urine

Distribution of Soil Transmitted Helminths


≈ 1.4 billion people infected




Trends in Parasitology

Figure 1. Global Distribution of Soil-Transmitted Helminths, 2010. Data from the *Global Atlas of Helminth Infection* were sourced to derive global estimates of soil-transmitted helminths (*Ascaris lumbricoides*, *Trichuris trichiura*, *Necator americanus*, and *Ancylostoma duodenale*) [83]. Reproduced, with permission, from [83].

PCR vs microscopy – what is the truth prevalence?




PCR vs microscopy – what is the truth prevalence?



Kaisar *et al* (submitted)

Community-based surveys



(* = not tested at each survey)

LUMC multiplex real-time PCR 2007-2016 (unpublished data)

Welke helminth targets in een feces PCR?

Klinische setting?

- *Strongyloides stercoralis*
- *Schistosoma* (*S. mansoni*)

Survey setting?

- *Ancylostoma duodenale*, *Necator americanus*
- *Ascaris lumbricoides*, *Trichuris trichiura*

Internationale markt voor een EQAS?

- Ja, ook in STH endemische gebieden
- Klinische materialen, maar ook DNA
- Feces in ethanol; stabiel voor opslag en transport

Pilot for helminth PCR EQAS – 2017/2018

HEMQAS

- Consortium of 18 academic, clinical, and public health laboratories and organizations from both endemic and non-endemic STH countries
- Clinical material from the field (mainly Gent)
 - Microscopy positive
 - Additional targets cannot be excluded
 - Negative material also included in the panel
- Central lab (NL) & 5 reference laboratories (NL 2x, USA 2x, Australia)
- Reproducibility, stability (12 feces in ethanol, 8 DNA)

Partner(s) contributing samples, microscopy confirmed, only human, mixed with ethanol

Early 2017

Spring 2017

Central Lab

- Homogenizing material
- Blinding; aliquotting
- Coordination of internal validation
- Selection of samples for reference labs
- ICT system, reporting

Internal validation

- Test feces samples & DNA
- Detectability
- Specificity
- Reproducibility feces (5x isolation)

Stage 1: Summer 2017

Reference partners

- 12 feces samples
- Sensitivity
- Specificity
- Reproducibility feces (5x isolation)
- Stability of target
- 4 species DNA (low, high)
- International transport issues
- Advisory role: selection of final panel

Stage 2: Spring 2018

Qbase Report


Central SKML office

- Qbase software
- reporting

- EQAS participant
- EQAS participant
- EQAS participant
- EQAS participant
- EQAS participant
- EQAS participant
- EQAS participant
- EQAS participant


- ➔ Blinded Samples
- ➔ Reporting via Excel
- ➔ Reporting via QBase

stool samples



DNA


DNA extract
worm/larvae
High concentration



- Reference laboratory 1
- Reference laboratory 2
- Reference laboratory 3
- Reference laboratory 4
- Reference laboratory 5

dilution

Low concentration



- Reference laboratory 1
- Reference laboratory 2
- Reference laboratory 3
- Reference laboratory 4
- Reference laboratory 5

Pilot for helminth PCR EQAS – 2017/2018

Final judgement per sample

- **Positive** and reproducible for target X:
 - All 5-plo positive and within Cq range ($X \pm 2SD$ or $2Ct$)
 - Positive by all 5 reference labs
- **Negative** for target Y:
 - All 5-plo negative
 - Negative by all 5 reference labs
- **Educational** for target Z:
 - Not uniformly positive (low loads)

Reference labs

Each of 12 stool and 8 DNA samples suitable to be send out

Clear what is **positive**, **negative**, **educational**

Illustrative examples: sample A

- Microscopy: species X + species Y
- PCR: **100% species X** + **60% species Y**
- PCR: **100% species V** + **100% species W**; microscopy missed
- PCR: **0% species S** + **0% species T**
- All 5-plo accordingly

Outcome of stage 2

Pilot for helminth PCR EQAS – 2017/2018

Participating labs (including reference labs)

Received 12 stool (ethanol) and 8 DNA samples

Extensive Qbase questionnaire about used procedures

Summer 2018: analysis of data



Image Credit: istockphoto.com/Feverpitched